

PRODUCT LINE



Corintha

Tuscan

Iona

Dora

Concentration

23 mg/ml

23 mg/ml

23 mg/ml

23 mg/ml

Cross-Linking



Gel Density

Super Soft

Soft

Medium-Soft

Medium-High

Particle Size

180-250 μ m

180-250 μ m

250-500 μ m

500-900 μ m

Lifting Capacity

Low

Moderate

Moderate

High

Indication

Superficial wrinkles

Lip contour

Moderate-to-deep wrinkles

Deep wrinkles

Perioral lines

Lip augmentation

Facial shaping
Nasolabial folds and
Marionette lines

Facial shaping

Periorbital lines

Superficial wrinkles

Lip augmentation

Cheeks
Deep nasolabial
folds and

Marionette lines

Injection Area

Middle Dermis

Middle dermis
Vermilion border

Middle Dermis

Deep Dermis

Package Content

2*1ml

2*1ml

2*1ml

2*1ml

Recommended Needle Type

30 G

30 G

27 G

27 G